University of Babylon College of Science Department of Chemistry Course No. Chsc 424



Undergraduate Studies
Physical chemistry
Fourth year - Semester 2
Credit Hour: 3 hrs.
Scholar units: three units

Lectures of Molecular Spectroscopy Second Semester, Scholar year 2024-2025 Prof. Dr. Abbas A-Ali Drea

\_\_\_\_\_

## **Lecture No. 8: Electron Spin Resonance Spectroscopy**

#### **Introduction:**

- A chemical species that contained one or more than one unpaired electron can be studied through Electron spin resonance technique.
- The radicals produced by such radiation damage, also can be studied to gives information about the locations of unpaired electrons and their mechanisms.
- ESR is very selective method of free radical identifications (neutral and negative radicals generally) due to the ability for detection at short life time of radicals (0.1 microsecond) and very low concentration approximately 10<sup>-10</sup> molar concentrations.

## Q// what the types and kinds of chemical species, ESR can detect?

ESR can detect the chemical species as follow:

- Transition metal ions and complexes Mn<sup>2+</sup>, Cu<sup>2+</sup>, Gd <sup>3+</sup> etc.
- Simple inorganic compounds: O<sub>2</sub>, NO, NO<sub>2</sub> ....
- Short-lived intermediate radicals \*OH, \*H, \*F etc. in kinetics study.
- Defects in crystals.
- Electrons trapped in radiation damaged sites.
- Stable organic radicals.
- Triplet states.

# Q// what are the Biological applications of ESR?

The biological application of ESR are-

- i- Paramagnetic cofactors: iron sulfur, copper proteins.
- ii- Free radicals of biological origin and their spin-trapping products.
- iii- Spin-labeling.

Therefore, this technique finds application in physics, chemistry, biology and medicine.

## 2- Physical concepts of ESR:

The basic physical concepts of EPR are analogous to those of nuclear magnetic resonance (NMR), since electron spins that are excited instead of spins of atomic nuclei.

When the molecules of a solid exhibit Para magnetism because of the unpaired electron spins, transitions can induced between spin states by applying a magnetic field and then supplying electromagnetic energy, usually in the microwave range of frequencies. The resulting absorption spectra are described as electron spin resonance (ESR). The interaction of an external magnetic field with an electron spin depends upon the magnetic moment associated with the spin, and the nature of an isolated electron spin is such that two and only two orientations are possible.

The application of the magnetic field then provides a magnetic potential energy that splits the spin states by an amount proportional to the magnetic field (Zeeman Effect), and then radio frequency radiation of the appropriate frequency can cause a transition from one spin state to the other. The energy associated with the transition is expressed in terms of the applied magnetic field B, the electron spin g-factor g, and the constant  $\mu$ . B which is called the Bohr magneton.

The difference between the energies of two possible spin states of a free electron in a magnetic field of flux density  $B_o$  is

$$\Delta E = \mu . \boldsymbol{B}. g. B_o - - - 1$$

Where  $\mu$ .B is Bohr's magnetron and g the Landé factor (almost exactly equal to 2 for a free electron) or Zeeman Effect. The absorption of electromagnetic radiation at the resonance frequency

$$\mathbf{v} = \Delta E/h$$
-----2

Equation 2 represent the Induces transitions from the low-energy level to the higher level. The numerical equation value is  $v = 2.8026 \cdot 10^{10} \cdot Bo$  (Hz / T). Due to the interaction between spin and lattice, the occupation of the higher level is again reduced. The absorbed energy is dissipated (radiated) to the environment in the form of heat. The resonance is determined by linearly varying the flux density B of the magnetic field over a small interval approximately the resonance value  $B_0$  and by recording the variation of the absorbed energy with Bas the signal of a so-called ESR spectrometer. Following figure represented the mechanism and signal of ESR.





Figure. Essentials of the EPR instrument.


## **Hyperfine interactions:**

The unpaired electron, which gives us the EPR spectrum, is very sensitive to local fields in its surroundings. Local fields arising from magnetic nuclei are permanent and independent of H. Interaction with neighboring nuclear

magnetic dipoles gives the nuclear hyperfine interaction and hyperfine splitting a corresponds to the NMR coupling constant J. A splitting are independent of the external field. For several equivalent nuclei n,  $(2n_MI_M + 1)$  transitions are observed for a nucleus M with a spin I. The relative intensities are given by Pascal's triangle for  $I = \frac{1}{2}$ . As in following rule:

Following figure, show this phenomenon. The phenomena depend also on the polarity of medium.






Hyperfine Splitting T Increases with Polarity of Medium, as example:


# Q/Explain the effect of medium polarity on the Hyperfine Splitting T Increases, In the ESR technique.

Answer/ as example for this phenomena in ESR, the nitroxide group is in equilibrium (resonance) between two structures, having the unpaired electron - on the O-atom (left, no hyperfine interaction) or N-atom (right, hyperfine interaction). A polar medium stabilizes the latter (electron is present on nitrogen atom) and thus increases the hyperfine splitting T.




**Applications:** 

1-Cyclooctatetraen anion: Observation of the spectrum shows that eight protons are equivalent1:8:28:56:70:56:28:8:1.



2-Butadien ion in liquid NH<sub>3</sub>: Two sets of equivalent protons: 2 and 4.



3-Pyrazine anion: two state with two different spectrum.

First with  $Na^+$  is the counter ion, second with  $K^+$  is the counter ion.

